

Sedens 410

Pressure Relieving Seat Cushion

Providing 24 Hour Pressure Area Care for Patient

- Alternating cushion works with Inspire dynamic mattress system pumps.
- Hybrid design with alternating therapy to relieve pressure from patient's buttock whilst reducing the risk of bottoming-out.

Effective seating solution to provide great pressure relief to patient when sitting in chair.

Crib 7 fire retardant, high moisture vapour transmission stretch cover with welded seams to wick moisture away from the surface while preventing any fluid ingress.



6 Air-cell-like one-piece pad with foam encased in to reduce risk of bottoming out.

Non-slip base prevents the cushion from sliding out of position, giving peace of mind to the user.

Specifications

Cushion	Cushion Dimension	44 x 41 x 9 cm
	Cell Height	3"
	Weight	1.11 kg
	Cover Material	Stretch PU Polyester
	Cell Material	Nylon PU
	Base material	Non-slip PU Polyester
	Fire Retardant Rating	EN 597-1 & EN 597-2
	Maximum Patient Weight	No weight limit, yet the cushion must be properly sized to the individual

Specifications are subject to change for improvement without notice

Inspire Community Trust

Unit 21, Optima Park,
Thames Road, Crayford,
Kent, DA1 4QX

T: 01322 520560
enquiry@inspiremedicare.org
www.inspirecommunitytrust.org

All rights reserved Specifications may change without notice.

Inspire Community Trust are a distributor of Apex Medical products.

CE